


Pipette MD Provides Fast, Reliable Calibrations For All Of Your Pipettes



Pipette MD offers:

- Trained professional technicians and support staff
- Calibration levels to meet your individual needs
- Calibration alerts—we will remind you of the next scheduled service interval
- Cost savings through a single service and calibration contact
- Service and Calibration in accordance with ISO 8655

Our Calibration Lab's skilled in-house technicians and support staff will collectively ensure fast, reliable service from receipt to delivery of your pipettes. Pipette MD service technicians are highly-trained professionals who are able to calibrate all makes and models of mechanical and electronic pipettes



Calibrates and maintains
all makes – all models – all volumes



“AS FOUND” CALIBRATION ONLY

Performance evaluation is conducted at two test volumes. Maximum and 10% of maximum with 3 measurements at each volume. This service includes “As Found” data only. All services are traceable to NIST and include a calibration report and label.

**Note: Multi-channel units are tested at center channel only. Visual linearity check is performed on all multi-channel units.*

Single-Channel	8-Channel	12-Channel
\$20.00	\$25.00	\$25.00

“AFTER SERVICE” CALIBRATION

Performance evaluation is conducted at two test volumes. 10% and 100% of maximum with 3 measurements at each volume. This service includes cleaning of the tip cone(s) and piston(s) with lubrication of pistons as well as “After Service” data only. All services are traceable to NIST and include a calibration report and label.

**Note: Multi-Channel units are tested at the center channel only. Visual linearity check is performed on all multi-channel units.*

Single-Channel	8-Channel	12-Channel
\$30.00	\$60.00	\$60.00

*Prices are available for a limited time and may change without notice. Please reference promotional code **CAL-11A** when ordering.*

Please Note: Only pipettes which are operational upon receipt can be calibrated. Non-operational pipettes, or pipettes which we cannot calibrate, will require additional labor and possibly spare parts. Customer will be notified via a quotation prior to additional work being performed.

Contact Us:

Pipette MD Technical Services

Web: www.pipetتمد.com

Email: info@pipetتمد.com

Phone: 1-877-MD-PIPET

Fax: 1-732-922-0557



Calibrates and maintains
all makes – all models – all volumes

