

BIOHIT

Limited Time Offer!

Biohit/Pipette MD 4-Step Preventative Maintenance Programs minimize the risk of pipetting errors and prolong the lifespan of your pipette



4-Step Mechanical Pipette Program

- Exterior cleaning of pipette
- Replacement of pistons and tip cones
- Leak and Linearity test
- Calibration to manufacturer specifications

PARTS & CALIBRATION AS LOW AS \$62.00

Preventative Maintenance Programs are available at highly discounted prices for all Biohit mechanical and electronic pipettes

PARTS & CALIBRATION AS LOW AS \$99.00

4-Step Electronic Pipette Program

- Exterior cleaning of pipette
- Replacement of battery, pistons, and tip cones
- Leak and Linearity test
- Calibration to manufacturer specifications



Once-a-year replacement of critical parts (tip cones, pistons, and battery for electronic models) with a basic calibration keeps your pipette ahead of the wear and tear curve



PREVENTATIVE MAINTENANCE PROGRAM

Preventative Maintenance Program

Programs are available for all Biohit brand pipettes

All pipettes are precision instruments, with mechanical and/or electronic components that are naturally subject to wear and possibly even damage over time, for example:

- The piston surface can deteriorate, wear, or corrode over time from friction and/or chemical attack.
- Tip cones are subject to chemical corrosion and wear, especially in the tip sealing area.
- Batteries within electronic pipettes have a normal lifespan of 1-2 years and must be replaced periodically.



Pipette failure rates can be significantly reduced with the application of a regular Preventative Maintenance Program. Enhance the reliability of your pipette and the integrity of your data!

Sartorius Biohit & Pipette MD recommend Preventative Maintenance programs be applied once a year.

Pipette Type	PM Program Includes	PM Program Cost	Cost if Purchased Separately
1-channel mLINE	Tip cone, piston, calibration	\$62.00	Over a \$130.00 value
8-channel mLINE	Tip cones, pistons, calibration	\$121.00	Over a \$400.00 value
12-channel mLINE	Tip cones, pistons, calibration	\$142.00	Over a \$600.00 value
1-channel eLINE	Tip cone, piston, battery, calibration	\$99.00	Over a \$190.00 value
8-channel eLINE	Tip cones, pistons, battery, calibration	\$157.00	Over a \$575.00 value
12-channel eLINE	Tip cones, pistons, battery, calibration	\$189.00	Over a \$760.00 value

Sartorius Biohit/Pipette MD can provide Preventative Maintenance Programs to any Biohit pipette, including Proline Plus, Proline Electronic, and ePET models. Contact us for additional details.



For more information:

Biohit Technical Services/Pipette MD

www.biohit.com/us

info@pipettemd.com

1-877-637-4738

- Offer expires December 31, 2012
- Pricing subject to change without notice
- May not be combined with any other special offers or promotions
- Parts-only pricing without calibration can be quoted separately at a non-promotional level upon request

